

Current CANDOR Catalog Number	Current CANDOR Product Name	Current CANDOR Pack Size-Specific Article Number	Pack Size	New Catalog Number	New Product Name	New Pack Size-Specific Article Number	New Product Description Including Pack Size Information
100	LowCross-Buffer OEM	100902-02	125 mL Sample	CR100	LowCross-Buffer	CR100-02	LowCross-Buffer 125 mL
100	LowCross-Buffer OEM	100902-04	1 L	CR100	LowCross-Buffer	CR100-04	LowCross-Buffer 1 L
100	LowCross-Buffer OEM	100902-06	10 L	CR100	LowCross-Buffer	CR100-06	LowCross-Buffer 10 L
101	LowCross-Buffer MILD OEM	101902-02	125 mL Sample	CR101	LowCross-Buffer MILD	CR101-02	LowCross-Buffer MILD 125 mL
101	LowCross-Buffer MILD OEM	101902-04	1 L	CR101	LowCross-Buffer MILD	CR101-04	LowCross-Buffer MILD 1 L
101	LowCross-Buffer MILD OEM	101902-06	10 L	CR101	LowCross-Buffer MILD	CR101-06	LowCross-Buffer MILD 10 L
102	LowCross-Buffer STRONG OEM	102902-02	125 mL Sample	CR102	LowCross-Buffer STRONG	CR102-02	LowCross-Buffer STRONG 125 mL
102	LowCross-Buffer STRONG OEM	102902-04	1 L	CR102	LowCross-Buffer STRONG	CR102-04	LowCross-Buffer STRONG 1 L
102	LowCross-Buffer STRONG OEM	102902-06	10 L	CR102	LowCross-Buffer STRONG	CR102-06	LowCross-Buffer STRONG 10 L
105	Sample Buffer OEM	105902-02	125 mL Sample	CR105	Sample Buffer	CR105-02	Sample Buffer 125 mL
105	Sample Buffer OEM	105902-04	1 L	CR105	Sample Buffer	CR105-04	Sample Buffer 1 L
105	Sample Buffer OEM	105902-06	10 L	CR105	Sample Buffer	CR105-06	Sample Buffer 10 L
106	Sample Buffer for IHC OEM	106902-02	125 mL Sample	CR106	Sample Buffer for IHC	CR106-02	Sample Buffer for IHC 125 mL
106	Sample Buffer for IHC OEM	106902-04	1 L	CR106	Sample Buffer for IHC	CR106-04	Sample Buffer for IHC 1 L
106	Sample Buffer for IHC OEM	106902-06	10 L	CR106	Sample Buffer for IHC	CR106-06	Sample Buffer for IHC 10 L
110	The Blocking Solution OEM	110902-02	125 mL Sample	CR110	The Blocking Solution	CR110-02	The Blocking Solution 125 mL
110	The Blocking Solution OEM	110902-04	1 L	CR110	The Blocking Solution	CR110-04	The Blocking Solution 1 L
110	The Blocking Solution OEM	110902-06	10 L	CR110	The Blocking Solution	CR110-06	The Blocking Solution 10 L
112	PlateBlock OEM	112902-02	125 mL Sample	CR112	PlateBlock	CR112-02	PlateBlock 125 mL
112	PlateBlock OEM	112902-04	1 L	CR112	PlateBlock	CR112-04	PlateBlock 1 L
112	PlateBlock OEM	112902-06	10 L	CR112	PlateBlock	CR112-06	PlateBlock 10 L
113	SmartBlock OEM	113902-02	125 mL Sample	CR113	SmartBlock	CR113-02	SmartBlock 125 mL
113	SmartBlock OEM	113902-04	1 L	CR113	SmartBlock	CR113-04	SmartBlock 1 L
113	SmartBlock OEM	113902-06	10 L	CR113	SmartBlock	CR113-06	SmartBlock 10 L
120	Coating Buffer pH 7.4 OEM	120902-02	125 mL Sample	CR120	Coating Buffer pH 7.4	CR120-02	Coating Buffer pH 7.4 125 mL
120	Coating Buffer pH 7.4 OEM	120902-04	1 L	CR120	Coating Buffer pH 7.4	CR120-04	Coating Buffer pH 7.4 1 L
120	Coating Buffer pH 7.4 OEM	120902-06	10 L	CR120	Coating Buffer pH 7.4	CR120-06	Coating Buffer pH 7.4 10 L
121	Coating Buffer pH 9.6 OEM	121902-02	125 mL Sample	CR121	Coating Buffer pH 9.6	CR121-02	Coating Buffer pH 9.6 125 mL
121	Coating Buffer pH 9.6 OEM	121902-04	1 L	CR121	Coating Buffer pH 9.6	CR121-04	Coating Buffer pH 9.6 1 L
121	Coating Buffer pH 9.6 OEM	121902-06	10 L	CR121	Coating Buffer pH 9.6	CR121-06	Coating Buffer pH 9.6 10 L
130	Antibody Stabilizer TRIS OEM	130902-02	125 mL Sample	CR130	Antibody Stabilizer TRIS	CR130-02	Antibody Stabilizer TRIS 125 mL
130	Antibody Stabilizer TRIS OEM	130902-04	1 L	CR130	Antibody Stabilizer TRIS	CR130-04	Antibody Stabilizer TRIS 1 L
130	Antibody Stabilizer TRIS OEM	130902-06	10 L	CR130	Antibody Stabilizer TRIS	CR130-06	Antibody Stabilizer TRIS 10 L
131	Antibody Stabilizer PBS OEM	131902-02	125 mL Sample	CR131	Antibody Stabilizer PBS	CR131-02	Antibody Stabilizer PBS 125 mL
131	Antibody Stabilizer PBS OEM	131902-04	1 L	CR131	Antibody Stabilizer PBS	CR131-04	Antibody Stabilizer PBS 1 L
131	Antibody Stabilizer PBS OEM	131902-06	10 L	CR131	Antibody Stabilizer PBS	CR131-06	Antibody Stabilizer PBS 10 L
140	Washing Buffer PBS 10x OEM	140902-02	125 mL Sample	CR140	Washing Buffer PBS 10x	CR140-02	Washing Buffer PBS 10x 125 mL
140	Washing Buffer PBS 10x OEM	140902-04	1 L	CR140	Washing Buffer PBS 10x	CR140-04	Washing Buffer PBS 10x 1 L
140	Washing Buffer PBS 10x OEM	140902-06	10 L	CR140	Washing Buffer PBS 10x	CR140-06	Washing Buffer PBS 10x 10 L
141	Washing Buffer PBS 10x for IHC OEM	141902-02	125 mL Sample	CR141	Washing Buffer PBS 10x for IHC	CR141-02	Washing Buffer PBS 10x for IHC 125 mL
141	Washing Buffer PBS 10x for IHC OEM	141902-04	1 L	CR141	Washing Buffer PBS 10x for IHC	CR141-04	Washing Buffer PBS 10x for IHC 1 L
141	Washing Buffer PBS 10x for IHC OEM	141902-06	10 L	CR141	Washing Buffer PBS 10x for IHC	CR141-06	Washing Buffer PBS 10x for IHC 10 L
145	Washing Buffer TRIS 10x OEM	145902-02	125 mL Sample	CR145	Washing Buffer TRIS 10x	CR145-02	Washing Buffer TRIS 10x 125 mL
145	Washing Buffer TRIS 10x OEM	145902-04	1 L	CR145	Washing Buffer TRIS 10x	CR145-04	Washing Buffer TRIS 10x 1 L
145	Washing Buffer TRIS 10x OEM	145902-06	10 L	CR145	Washing Buffer TRIS 10x	CR145-06	Washing Buffer TRIS 10x 10 L
146	Washing Buffer TRIS 10x for IHC OEM	146902-02	125 mL Sample	CR146	Washing Buffer TRIS 10x for IHC	CR146-02	Washing Buffer TRIS 10x for IHC 125 mL
146	Washing Buffer TRIS 10x for IHC OEM	146902-04	1 L	CR146	Washing Buffer TRIS 10x for IHC	CR146-04	Washing Buffer TRIS 10x for IHC 1 L
146	Washing Buffer TRIS 10x for IHC OEM	146902-06	10 L	CR146	Washing Buffer TRIS 10x for IHC	CR146-06	Washing Buffer TRIS 10x for IHC 10 L
160	Liquid Plate Sealer OEM	160902-02	125 mL Sample	CR160	Liquid Plate Sealer	CR160-02	Liquid Plate Sealer 125 mL
160	Liquid Plate Sealer OEM	160902-04	1 L	CR160	Liquid Plate Sealer	CR160-04	Liquid Plate Sealer 1 L
160	Liquid Plate Sealer OEM	160902-06	10 L	CR160	Liquid Plate Sealer	CR160-06	Liquid Plate Sealer 10 L
163	Liquid Plate Sealer animal-free OEM	163902-02	125 mL Sample	CR163	Liquid Plate Sealer animal-free	CR163-02	Liquid Plate Sealer animal-free 125 mL
163	Liquid Plate Sealer animal-free OEM	163902-04	1 L	CR163	Liquid Plate Sealer animal-free	CR163-04	Liquid Plate Sealer animal-free 1 L
163	Liquid Plate Sealer animal-free OEM	163902-06	10 L	CR163	Liquid Plate Sealer animal-free	CR163-06	Liquid Plate Sealer animal-free 10 L
180	Assay Defender OEM	180902-02	125 mL Sample	CR180	Assay Defender	CR180-02	Assay Defender 125 mL
180	Assay Defender OEM	180902-04	1 L	CR180	Assay Defender	CR180-04	Assay Defender 1 L
180	Assay Defender OEM	180902-06	10 L	CR180	Assay Defender	CR180-06	Assay Defender 10 L
222	HRP-Protector OEM	222902-02	125 mL Sample	CR222	HRP-Protector	CR222-02	HRP-Protector 125 mL
222	HRP-Protector OEM	222902-04	1 L	CR222	HRP-Protector	CR222-04	HRP-Protector 1 L
222	HRP-Protector OEM	222902-06	10 L	CR222	HRP-Protector	CR222-06	HRP-Protector 10 L
270	LowCross HRP-Stab OEM	270902-02	125 mL Sample	CR270	LowCross HRP-Stab	CR270-02	LowCross HRP-Stab 125 mL
270	LowCross HRP-Stab OEM	270902-04	1 L	CR270	LowCross HRP-Stab	CR270-04	LowCross HRP-Stab 1 L
270	LowCross HRP-Stab OEM	270902-06	10 L	CR270	LowCross HRP-Stab	CR270-06	LowCross HRP-Stab 10 L
271	LowCross HRP-Stab MILD OEM	271902-02	125 mL Sample	CR271	LowCross HRP-Stab MILD	CR271-02	LowCross HRP-Stab MILD 125 mL
271	LowCross HRP-Stab MILD OEM	271902-04	1 L	CR271	LowCross HRP-Stab MILD	CR271-04	LowCross HRP-Stab MILD 1 L
271	LowCross HRP-Stab MILD OEM	271902-06	10 L	CR271	LowCross HRP-Stab MILD	CR271-06	LowCross HRP-Stab MILD 10 L
272	LowCross HRP-Stab STRONG OEM	272902-02	125 mL Sample	CR272	LowCross HRP-Stab STRONG	CR272-02	LowCross HRP-Stab STRONG 125 mL
272	LowCross HRP-Stab STRONG OEM	272902-04	1 L	CR272	LowCross HRP-Stab STRONG	CR272-04	LowCross HRP-Stab STRONG 1 L
272	LowCross HRP-Stab STRONG OEM	272902-06	10 L	CR272	LowCross HRP-Stab STRONG	CR272-06	LowCross HRP-Stab STRONG 10 L
300	LowCross-Buffer protein-free OEM	300902-02	125 mL Sample	CR300	LowCross-Buffer protein-free	CR300-02	LowCross-Buffer protein-free 125 mL
300	LowCross-Buffer protein-free OEM	300902-04	1 L	CR300	LowCross-Buffer protein-free	CR300-04	LowCross-Buffer protein-free 1 L
300	LowCross-Buffer protein-free OEM	300902-06	10 L	CR300	LowCross-Buffer protein-free	CR300-06	LowCross-Buffer protein-free 10 L
301	LowCross-Buffer MILD protein-free OEM	301902-02	125 mL Sample	CR301	LowCross-Buffer MILD protein-free	CR301-02	LowCross-Buffer MILD protein-free 125 mL
301	LowCross-Buffer MILD protein-free OEM	301902-04	1 L	CR301	LowCross-Buffer MILD protein-free	CR301-04	LowCross-Buffer MILD protein-free 1 L
301	LowCross-Buffer MILD protein-free OEM	301902-06	10 L	CR301	LowCross-Buffer MILD protein-free	CR301-06	LowCross-Buffer MILD protein-free 10 L
302	LowCross-Buffer STRONG protein-free OEM	302902-02	125 mL Sample	CR302	LowCross-Buffer STRONG protein-f	CR302-02	LowCross-Buffer STRONG protein-free 125 mL
302	LowCross-Buffer STRONG protein-free OEM	302902-04	1 L	CR302	LowCross-Buffer STRONG protein-f	CR302-04	LowCross-Buffer STRONG protein-free 1 L
302	LowCross-Buffer STRONG protein-free OEM	302902-06	10 L	CR302	LowCross-Buffer STRONG protein-f	CR302-06	LowCross-Buffer STRONG protein-free 10 L
305	Sample Buffer protein-free OEM	305902-02	125 mL Sample	CR305	Sample Buffer protein-free	CR305-02	Sample Buffer protein-free 125 mL
305	Sample Buffer protein-free OEM	305902-04	1 L	CR305	Sample Buffer protein-free	CR305-04	Sample Buffer protein-free 1 L
305	Sample Buffer protein-free OEM	305902-06	10 L	CR305	Sample Buffer protein-free	CR305-06	Sample Buffer protein-free 10 L
322	HRP-Protector protein-free OEM	322902-02	125 mL Sample	CR322	HRP-Protector protein-free	CR322-02	HRP-Protector protein-free 125 mL
322	HRP-Protector protein-free OEM	322902-04	1 L	CR322	HRP-Protector protein-free	CR322-04	HRP-Protector protein-free 1 L
322	HRP-Protector protein-free OEM	322902-06	10 L	CR322	HRP-Protector protein-free	CR322-06	HRP-Protector protein-free 10 L
330	Antibody Stabilizer TRIS protein-free OEM	330902-02	125 mL Sample	CR330	Antibody Stabilizer TRIS protein-fre	CR330-02	Antibody Stabilizer TRIS protein-free 125 mL
330	Antibody Stabilizer TRIS protein-free OEM	330902-04	1 L	CR330	Antibody Stabilizer TRIS protein-fre	CR330-04	Antibody Stabilizer TRIS protein-free 1 L
330	Antibody Stabilizer TRIS protein-free OEM	330902-06	10 L	CR330	Antibody Stabilizer TRIS protein-fre	CR330-06	Antibody Stabilizer TRIS protein-free 10 L
331	Antibody Stabilizer PBS protein-free OEM	331902-02	125 mL Sample	CR331	Antibody Stabilizer PBS protein-free	CR331-02	Antibody Stabilizer PBS protein-free 125 mL
331	Antibody Stabilizer PBS protein-free OEM	331902-04	1 L	CR331	Antibody Stabilizer PBS protein-free	CR331-04	Antibody Stabilizer PBS protein-free 1 L
331	Antibody Stabilizer PBS protein-free OEM	331902-06	10 L	CR331	Antibody Stabilizer PBS protein-free	CR331-06	Antibody Stabilizer PBS protein-free 10 L